


**Product Name: Recombinant Human ZUP1 Protein**  
**Catalog #: PEH10024**



## Summary

<b>Name</b>	ZUP1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96AP4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	31 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	AHFEQNTLER NFERINTVQY GTSDNKKDNT LQCGMEVNSS ILSGCASNHP KNSAQNLTKD STLKHEGFYS ENLTESRKFL KSREKQSSLT EIKGSVYETT YSPPECPFCG KIEEHSEDEME THVKTKHANL LDIPLEDQDQ PLYDCPMCGL ICTNYHILQE HVDLHLEENS FQQGMDRVQC SGDLQLAHQL QQEEDRKRRS EESRQEIEEF QKLQRQYGLD NSGGYKQQQL RNMEIEVNRG RMPPSEFHRR

## image



## Background

<b>Alternative Names</b>	Zinc finger-containing ubiquitin peptidase 1 (EC 3.4.19.12) (Lys-63-specific deubiquitinase ZUFSP) (DUB) (Zinc finger with UFM1-specific peptidase domain protein)
--------------------------	--

**Product Name: Recombinant Human ZUP1 Protein**  
**Catalog #: PEH10024**



---

**Note**

For Research Use Only , Not for Diagnostic Use