


Product Name: Recombinant Human ZSCAN21 Protein
Catalog #: PEH10027



Summary

Name	ZSC21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5A6
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEEKGKYLPS LEMFRQRFRQ FGYHDTGPR EALSQRVLC CEWLRPEIHT KEQILELLVL EQFLTLPQE LQAWVQEHCP ESAEEAVTLL EDLERELDEP GHQVSTPPNE QKPVWEKISS SGTAKESPSS MQPQPLETSH KYESWGPLYI QESGEEQEFA QDPRKVRDCR LSTQHEESAD EQKGSEAEGL KGDIIISVIA NKPEASLERQ CVNLENEKGT

image



Background

Alternative Names	Zinc finger and SCAN domain-containing protein 21 (Renal carcinoma antigen NY-REN-21) (Zinc finger protein 38 homolog) (Zfp-38)
--------------------------	---

Note

Product Name: Recombinant Human ZSCAN21 Protein
Catalog #: PEH10027



For Research Use Only , Not for Diagnostic Use