


Product Name: Recombinant Human SP3 Protein
Catalog #: PEH10028



Summary

Name	SP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02447
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TQAQNIQVST AQPVVQHLQL QESQQPTSQA QIVQGITPQT IHGVQASGQN ISQQALQNLQ LQLNPGTFLI QAQTVTPSGQ VTWQTFQVQG VQNLQNLQIQ NTAAQQITLT PVQTLTLGQV AAGGAFSTP VSLSTGQLPN LQTVTVNSID SAGIQLHPGE NADSPADIRI KEEEPDPEEW QLSGDSTLNT NDLTHLRVQV VDEEGDQQHQ

image



Background

Alternative Names Transcription factor Sp3 (SPR-2)

Note

Product Name: Recombinant Human SP3 Protein
Catalog #: PEH10028



For Research Use Only , Not for Diagnostic Use