


Product Name: Recombinant Human SPATA5 Protein
Catalog #: PEH10038



Summary

Name	AFG2H_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NB90
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDMEINEEEL TGCILRKLDG KIVLPGNFLY CTFYGRPYKL QVLRVKGADG MILGGPQSDS DTDAQRMAFE QSSMETSSLE LSLQLQLDL EDTQIPTSR TPYKPIDDRI TNKASDVLLD VTQSPGDGSG LMLEEV TGLK CNFESAREGN EQLTEERLL KFSIGAKCNT DTFYFISSTT RVNFTEIDKN SKEQDNQFKV

image



Background

Alternative Names ATPase family protein 2 homolog (EC 3.6.4.10) (Spermatogenesis-associated factor protein) (Spermatogenesis-associated protein 5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SPATA5 Protein
Catalog #: PEH10038

