

Product Name: Recombinant Human EID1 Protein
Catalog #: PEH10041



Summary

Name	EID1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6B2
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YEESDLQMD VMPGEGDLPQ MEVGSGSREL SLRPSRSGAQ QLEEEGPMEE EEAQPMAAPE GKRLANGPN AGEQPGQVAG ADFESEDEGE EFDDWEDDYD YPREEQLSGA GYRVSAALEE ADKMFLRTRE PALDGGFQMH YEKTPFDQLA FIEELFSLMV VNRLTEELGC

image



Background

Alternative Names	EP300-interacting inhibitor of differentiation 1 (21 kDa pRb-associated protein) (CREBBP/EP300 inhibitory protein 1) (E1A-like inhibitor of differentiation 1) (EID-1)
--------------------------	--

Note

Product Name: Recombinant Human EID1 Protein
Catalog #: PEH10041



For Research Use Only , Not for Diagnostic Use