


Product Name: Recombinant Human RPL3 Protein
Catalog #: PEH10043



Summary

Name	RL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P39023
Host	E.coli
Species	Human
Molecular Mass	88 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGSLGFLPRK RSSRHRGKVK SFPKDDPSKP VHLTAFLGYK AGMTHIVREV DRPGSKVNKK EVVEAVTIVE TPPMVVVGIV GYVETPRGLR TFKTVFAEHI SDECKRRFYK NWHKSKKAF TKYCKKWQDE DGKKQLEKDF SSMKKYCQVI RVIAHTQMRL LPLRQKKAHL MEIQVNGGTV AEKLDWARER LEQQVPVNQV FGQDEMIDVI GVTKGKGYKG VTSRWHTKKL PRKTHRGLRK VACIGAWHPA RVAFSVARAG QKGYHHRTEI NKKIYKIGQG YLIKDGKLIK NNASTDYDLS DKSINPLGGF VHYGEVTNDF VMLKGCVVGT KKRVLTLRKS LLVQTKRRAL EKIDLKFIDT TSKFGHGRFQ TMEEKKAFMG

image



Background

Alternative Names 60S ribosomal protein L3 (HIV-1 TAR RNA-binding protein B) (TARBP-B)

Product Name: Recombinant Human RPL3 Protein
Catalog #: PEH10043



(Large ribosomal subunit protein uL3)

Note

For Research Use Only , Not for Diagnostic Use