


Product Name: Recombinant Human FBXW8 Protein
Catalog #: PEH10044



Summary

Name	FBXW8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N3Y1
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QAPKKRRRPE AAERRARRPE VGSGRGEQAS GDPALAQRLLEGAGRPPAAR ATRAEGQDVA SRSRSPLARE GAGGGEQLVD

image



Background

Alternative Names F-box/WD repeat-containing protein 8 (F-box and WD-40 domain-containing protein 8) (F-box only protein 29)

Note

For Research Use Only , Not for Diagnostic Use