


Product Name: Recombinant Human FMNL3 Protein
Catalog #: PEH10048



Summary

Name	FMNL3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVF7
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDSEALARVG PAELSEGMPP SDLDLLAPAP PEEVLPLPP PPAPPLPPPP PPLPDKCPPA PPLPGAAPSV VLTVGLSAIR IKKPIKTKFR LPVFNWTALK PNQISGTVFS ELDDEKILED LDLDKFEELF KTKAQGPALD LICSKNKTAQ KAASKVTLL E ANRAKNLAIT LRKAGRSAEE ICRAIHTFDL QTLPVDFVEC LMRFLPTEAE VKLLRQYERE RQPLEELAAE DRFMLLSKV ERLTQRMAGM AFLGNFQDNL QMLTPQLNAI IAASASVKSS QKLKQMLEII LALGNYMNSS KRGAVYGFKL QSLDLLLDTK STDRKMTLLH FIALTVKEY PDLANFWHEL HFVEKAAAVS LENVLLDVKE LGRGMELIRR ECSIHDNSVL RNFLSTNEGK LDKLQRDAKT AEEAYNAVVR YFGESPKTTP PSVFFPVFVR FIRSYKEAEQ ENEARKKQEE VMREKQLAQE AKKLDKTPS QRNKWQQQEL

image



Background

Product Name: Recombinant Human FMNL3 Protein
Catalog #: PEH10048



Alternative Names

Formin-like protein 3 (Formin homology 2 domain-containing protein 3)
(WW domain-binding protein 3) (WBP-3)

Note

For Research Use Only , Not for Diagnostic Use