


Product Name: Recombinant Human MPP4 Protein
Catalog #: PEH10050



Summary

Name	MPP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JB8
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PDKKDMKLST ATNPQNGLSQ ILRLVLQELS LFYGRDVNGV CLLYDLLHSP WLQALLKIYD CLQEFKEKKL VPATPHAQVL SYEVVELLRE TPTSPEIQEL RQMLQAPHFK ALLSAHDTIA QKDFEPLLPP LPDNIPESSE AMRIVCLVKN QQPLGATIKR HEMTGDILVA RIIHGGLAER SGLLYAGDKL VEVNGVSVEG LDPEQVIHIL AMSRGTIMFK VVPVSDPPVN SQQMYYVVRAM TEYWPQEDPD IPCMDAGLPF QKGDILQIVD QNDALWWQAR KISDPATCAG LVPSNHLLKR KQREFWWSQP YQPHTCLKST LSISMEEDD MKIDEKCVEA DEETFESEEL SEDKEEFVGY GQKFFIAGFR

image



Background

Alternative Names MAGUK p55 subfamily member 4 (Amyotrophic lateral sclerosis 2)

Product Name: Recombinant Human MPP4 Protein
Catalog #: PEH10050



chromosomal region candidate gene 5 protein) (Discs large homolog 6)

Note

For Research Use Only , Not for Diagnostic Use