


Product Name: Recombinant Human MMAA Protein
Catalog #: PEH10056



Summary

Name	MMAA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVH4
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKDHTEGLSD KEQRFVDKLY TGLIQGQRAC LAEAITLVES THSRKKELAQ VLLQKVLLYH REQEQSNK GK PLAFRVGLSG PPGAGKSTFI EYFGKMLTER GHKLSVLAVD PSSCTSGGSL LGDKTRMTEL SRDMNAYIRP SPTRGTLGGV TRTTNEAILL CEGAGYDIIL IETVGVGQSE FAVADMVDMF VLLPPAGGD ELQGIKRGII EMADLVAVTK

image



Background

Alternative Names Methylmalonic aciduria type A protein, mitochondrial (EC 3.6.-.-)

Note

Product Name: Recombinant Human MMAA Protein
Catalog #: PEH10056



For Research Use Only , Not for Diagnostic Use