


Product Name: Recombinant Human ANKRD55 Protein
Catalog #: PEH10058



Summary

Name	ANR55_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3KP44
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HQKDPSRDRY REEDTSEVND IITTFDSIVG TNCQEQPGDQ VAMVEFKKKT SDNSKYLLPE KKPLARKGLP PIRTQSLPPI TLGNNFLTAS HRATSHAGLS SAPHHMAQRS QKSRSEQDLL NNRTGCQMLL DNPWKSDSNQ VFSYKVWTVS

image



Background

Alternative Names Ankyrin repeat domain-containing protein 55

Note

For Research Use Only , Not for Diagnostic Use