


Product Name: Recombinant Human MAGEB16 Protein
Catalog #: PEH10061



Summary

Name	MAGBG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A2A368
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETQSLEVAQV SKALEKTLLS SSHPLVPGKL KEAPAAKAES PLEVQSFCS SSIAVTTTSS SESDEASSNQ EEEDSPSSSE DTSDPRNVPA DALDQKVAFL VNFMLHKCQM KKPITKADML KIIKDDESH FSEILLRASE HLEMIFGLDV VEVDPTTHCY GLFIKLGLTY DGMLSGEKGV PKTGLLIIVL GVIFMKGNRA TEEEVWEVLN LTGVYSGKKH FIFGEPRLI TKDFVKEKYL EYQQVANS DP ARYEFLWGPR AKAETSKMKV LEFVAKVHGS YPHSFPSQYA EALKEEEERA

image



Background

Alternative Names Melanoma-associated antigen B16 (MAGE-B16 antigen)

Product Name: Recombinant Human MAGEB16 Protein
Catalog #: PEH10061



Note

For Research Use Only , Not for Diagnostic Use