


Product Name: Recombinant Human AHI1-DT Protein
Catalog #: PEH10062



Summary

Name	CF217_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	POC7V0
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VKEAQKYGED SQKSHSPGTP GPRSSVTTLS ASALSDSSSP DTPRRGPGR PSTPARAPAT SAPMMYSRRG VRRTPARPAGA DTRSSANQLP QPSGACANAD SAPPADV SAC LRRRSHGDC VPRSRRRPRP

image



Background

Alternative Names Putative uncharacterized protein encoded by LINC00271

Note

For Research Use Only , Not for Diagnostic Use