


Product Name: Recombinant Human MAGEB6 Protein
Catalog #: PEH10065



Summary

Name	MAGB6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N7X4
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TCEKRQETNG QPQGLTGPQA TAEKQEESH SSSSRACLG DCCRSSDASI PQESQGVSP GSPDAVVSYS KSDVAANGQD EKSPSTRDA SVPQESQGAS PTGSPDAGVS

image



Background

Alternative Names Melanoma-associated antigen B6 (Cancer/testis antigen 3.4) (CT3.4) (MAGE-B6 antigen)

Note

For Research Use Only , Not for Diagnostic Use