


Product Name: Recombinant Human ZNG1A Protein
Catalog #: PEH10066



Summary

Name	CBWD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BRT8
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEDPAEEDCP ELVPIETTQS EEEEKSGLGA KIPVTIITGY LGAGKTLLN YILTEQHSCR VAVILNEFGE GSALEKSLAV SQGGELYEEW LELRNGCLCC SVKDSGLRAI ENLMQKKGKF

image



Background

Alternative Names	COBW domain-containing protein 1 (Cobalamin synthase W domain-containing protein 1) (COBP) (NPC-A-6 COBW domain-containing protein 1) (NPC-A-6)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNG1A Protein
Catalog #: PEH10066

