


Summary

Name	CFA58_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T655
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEESAFEEME RDFQGVLEL SGDKSLEKFR IEYERLHAVM KKSVDNEKRL MAKCRELNAE IVNSAKVAT ALKLSQDDQT TIASLKKEIE KAWKMOVDSAY DKEQKAKETI LALKEEIVNL TKLVEQGSGL SMDQHSNIRD LLRFKEEVTK ERDQLLSEVV KLRESLAQTT EQQQETERSK EEAHHAISQF QQEIQQRQNE ASREFRKKEK LEKELKQIQA DMDSRQTEIK ALQQYVQKSK EELQKLEQQL KEQKILNERA AKELEQFQMR NAKLQQENEQ HSLVCEQLSQ ENQQKALELK AKEEEVHQMR LDIGKLNKIR EQIHKKLHHT EDQKAEVEQH KETLKNQIVG LEREVEASKK QAELDRKAMD ELLRERDILN KNMLKAVNAT QKQTDLVKLN EQAKRNLEGE IQNYKDEAQQ QRKIIFHLEK ERDRYINQAS DLTQKVLMMN

image



Background

Product Name: Recombinant Human CFAP58 Protein
Catalog #: PEH10068



Alternative Names

Cilia- and flagella-associated protein 58 (Coiled-coil domain-containing protein 147)

Note

For Research Use Only , Not for Diagnostic Use