

Product Name: Recombinant Human KLHL33 Protein
Catalog #: PEH10071



Summary

Name	KLH33_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NCF5
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DWEEMAPLSQ ARSLFSLVAL DGKLYALGGR HNDVALDSVE TYNPELNVWR PAPALPAPCF AHAAAILEGQ LYVSGGCGGT GQYLASLMHY DPKLEKPGTF LSPMGVPRAG HVMAALGGRL YVAGGLGETE DLLSFEAYEL RTDSWTHLAP LPSPHVGAAS AVLQGELLVL

image



Background

Alternative Names Kelch-like protein 33

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KLHL33 Protein
Catalog #: PEH10071

