


Product Name: Recombinant Human ARL10 Protein
Catalog #: PEH10075



Summary

Name	ARL10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N8L6
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YFGRGRERRW DRGEAWWGAE AARLPEWDEW DPEDEEDEEP ALEELEQREV LVLGLDGAGK STFLRVLSGK PPLEGHIPTW GFNSVRLPTK DFEVDLLEIG GSQNLRFYWK EFVSEVDVLV FVVD SADRLR LPWARQELHK LLDKDPDLPV VVVANKQDLS EAMSMGELQR ELGLQAIDNQ

image



Background

Alternative Names ADP-ribosylation factor-like protein 10

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARL10 Protein
Catalog #: PEH10075

