


Product Name: Recombinant Human AGAP4 Protein
Catalog #: PEH10077



Summary

Name	AGAP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96P64
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSIPDEQLHS FAVSTVHIMK KRNGGSLNN YSSIPSTPS TSQEDPQFSV PPTANTPTPV CKRSMRWSNL FTSEKGSDDP KERKAPENHA

image



Background

Alternative Names Arf-GAP with GTPase, ANK repeat and PH domain-containing protein 4 (AGAP-4) (Centaurin-gamma-like family member 1) (Centaurin-gamma-like family member 5)

Note

For Research Use Only , Not for Diagnostic Use