

Product Name: Recombinant Human CCDC168 Protein
Catalog #: PEH10084



Summary

Name	CC168_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NDH2
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLDAHLEKEI KTEKEMHQPI PFTETIESV VSPIMELSHA ENVKSTQKTQ TDCKCTADSE TPSPISGKSL IGDPLNQTRE SYIPSNGSDT REMGYCFAEE KTEIPKDLPA TSPETFNYCT PVLSCSKVMK KRVTFALTTS TAKPKCVNTK AVKPSISETV SVTSHRKKSE LDFKTKFKKI NQTKGLVPEC LNTLCSPMHS RLQREFCLPA

image



Background

Alternative Names Coiled-coil domain-containing protein 168

Note

Product Name: Recombinant Human CCDC168 Protein
Catalog #: PEH10084



For Research Use Only , Not for Diagnostic Use