


Product Name: Recombinant Human FAM131C Protein
Catalog #: PEH10089



Summary

Name	F131C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96AQ9
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRSRPGNYNV AALATSSLVG VVQSIKDHIT KPTAMARGRV AHLIEWKGWS AQPAGWELSP AEDEHYCCLP DELREARFAA GVAEQFAITE ATLSAWSSLD EELHPENSP QGIVQLQDLE SIYLQDSLPS GPSQDDSLQA

image



Background

Alternative Names Protein FAM131C

Note

For Research Use Only , Not for Diagnostic Use