


Product Name: Recombinant Human FAM222B Protein
Catalog #: PEH10092



Summary

Name	F222B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WU58
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSTLAHPQAAQ ALARQQALQH AQTLAHAPPQ TLQHPQGIPP PQALSHPQSL QQPQGLGHPQ PMAQTQGLVH PQALAHQGLQ HPHNPLLHGG RKMPDSDAPP

image



Background

Alternative Names Protein FAM222B

Note

For Research Use Only , Not for Diagnostic Use