


Product Name: Recombinant Human GSN-AS1 Protein
Catalog #: PEH10095



Summary

Name	GSAS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRJ2
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGKATLLKSH SKQPAQPGWE PGIRAPSPVP ASSLQDHSRL TSLSRTGKEQ RRTLSLIRKT SGTPTTESTVA TAAASTTEVP

image



Background

Alternative Names	Putative uncharacterized protein GSN-AS1 (GSN antisense RNA 1) (GSN antisense gene protein 1) (Protein MOST-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use