


Product Name: Recombinant Human TTC9C Protein
Catalog #: PEH10098



Summary

Name	TTC9C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5M4
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YKEEGNQRYR EGKYRDAVSR YHRALLQLRG LDPSLPSPLP NLGPQGPALT PEQENILHTT QTDCYNNLAA CLLQMEPVNY ERVREYSQKV LERQPDNAKA LYRAGVAFFH LQDYDQARHY LLAAVNRQPK DANVRRYLQL TQSELSSYHR

image



Background

Alternative Names Tetratricopeptide repeat protein 9C (TPR repeat protein 9C)

Note

For Research Use Only , Not for Diagnostic Use