



Product Name: Recombinant Human C10orf82 Protein
Catalog #: PEH10099



Summary

Name	CJ082_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WW14
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TFQEKTRYK EQLRELCCAV ATAPKLKPVN SEETVLQALH QYNLQYHPLI LECKYVKKPL QEPIPGWAG YLPRAKVTEF GCGTRYTVMA KNCYKDFLEI TERAKKAHLK PYEEIYGVSS TKTSAPSPKV LQHEELLPKY PDFSIPDGSC PALGRPLRED

image



Background

Alternative Names Uncharacterized protein C10orf82

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C10orf82 Protein
Catalog #: PEH10099

