


Product Name: Recombinant Human C8orf76 Protein
Catalog #: PEH10103



Summary

Name	CH076_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96K31
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKSFFPHSGK DCLLCFPETL PESSLFSVEA NSSNSQKNEK ALTNIQNCMA EKRETVLIET QLKACASFIR TRLLLQFTQP QQTSFALERN LRTQQEIEDK MKGFSFKEDT LLLIAEVMGE DIPEKIKDEV HPEVKCVGSV ALTALVTVSS EEFEDKWFRK IKDHFPCFEN

image



Background

Alternative Names Uncharacterized protein C8orf76

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C8orf76 Protein
Catalog #: PEH10103

