



Summary

Name	CI043_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TAL5
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PALRYPERLK KLHNLKTEGY RKQQQRQQQQ QQQQKKVKTP IKKQEAKKKA KSDPGIQSTS HKHPVTTVHD RLYGYRTLPG QNSDMKQQQQ MEKGTTSKQD STERPKMNYD DHADFHHSVK SPELYETEPT NKDISAPVDA VPEAQAARQK KISFNFSEIM ASTGWNSELK LLRILQDTDD EDEEDQSSGA

image



Background

Alternative Names Uncharacterized protein C9orf43

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C9orf43 Protein
Catalog #: PEH10104

