



Product Name: Recombinant Human C1orf131 Protein
Catalog #: PEH10107



Summary

Name	CA131_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NDD1
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ILAAAVPPSS LKNNREQVEV VEFHSNKKRK LTPDHNKNTK QANPSVLERD VDTQEFNLEK ARLEVHRFGI TGYGKGKERI LEQERAIMLG AKPPKKSQYVN YKVLQEIQIKE KKAKEEEEKR LAQETDIFKK KKRKGQEDRK SKKKSAPSIL

image



Background

Alternative Names Uncharacterized protein C1orf131

Note

For Research Use Only , Not for Diagnostic Use