


Product Name: Recombinant Human GOLGA8IP Protein
Catalog #: PEH10109



Summary

Name	GOG8I_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NC78
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEISQLKNTI KSLKQKKQV EHQLKEEKKA NIKKQKAKRV LEVQIQTINI QKEELNTDLY HMKRSLRYFE EKSKDLAVRL QHSLQRKGEL ESVLSDVMAT QKKKANQLSS PSKAGTEWKL EQSMREEALL KVQLTQLKES FQQLQLERHE YAEHLKGERA RWQQRMRKMS QEICTLKKEK QQDMRRVEKL ERSLSKLKNQ MAEPLPPEPP AVPSEVELQH LRKELERVAG ALQAQVKNNQ RISLLNRGQE ERIREQEERL RKQEERIQEQ HKSLQQLAKP QSVFEPPNNE NKSALQLEQQ VKELQEKLGE EHLEAASQQN QQLTAQLSLM ALPGEGHGGE HLDSEGEAAP QPMPSVPEDL ESREAMSSFM DHLEEKADLS

image



Background

Alternative Names	Putative golgin subfamily A member 8I (Golgin subfamily A member 8I)
--------------------------	--



Product Name: Recombinant Human GOLGA8IP Protein
Catalog #: PEH10109

pseudogene)

Note

For Research Use Only , Not for Diagnostic Use