


Product Name: Recombinant Human DCUN1D4 Protein
Catalog #: PEH10112



Summary

Name	DCNL4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92564
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RHDSMYRKYD STRIKTEEEA FSSKRCLEWF YEYAGTDDVV GPEGMEKFCE DIGVEPENVV MLVLAWK LDA QNMGYFTLQE WLKGMTSLQC DTTEKL RNTL DYLR SFLNDS TNFKLIYRYA FDFAREKDQR SLDINTAKCM LGLLLGKIWP LFPVFHQFLE QSKYKVINKD QWCNVLEFSR TINL DLSNYD EDGAWPVLLD

image



Background

Alternative Names DCN1-like protein 4 (DCNL4) (DCUN1 domain-containing protein 4) (Defective in cullin neddylation protein 1-like protein 4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DCUN1D4 Protein
Catalog #: PEH10112

