


Product Name: Recombinant Human ANKRD24 Protein
Catalog #: PEH10114



Summary

Name	ANR24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TF21
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LILHLLQEAA QRPSPPSALT EDDSGEASSQ NSMSSHGKQG APKKRKAPPP PASIPMPDDR DAYEEIVRLR QERGRLLQKI RGLQHKERR QQESPEASSL HILERQVQEL QQLLVERQEE KESLGREVES LQSRLSLEN ERENTSYDVT TLQDEEGELP DLPGAEVLLS RQLSPSAQEH LASLQEQVAV LTRQNQELME KVQILENFEK DETQMEVEAL AEVIPLALYD SLRAEFDQLR RQHAEALQAL RQQETREVPR

image



Background

Alternative Names Ankyrin repeat domain-containing protein 24

Product Name: Recombinant Human ANKRD24 Protein
Catalog #: PEH10114



Note

For Research Use Only , Not for Diagnostic Use