


Product Name: Recombinant Human AGAP9 Protein
Catalog #: PEH10115



Summary

Name	AGAP9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VTM2
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NSKPLGSEAI DNQGVSLFED QQQGSVCPSE SEIYEAGAED RMAGAPMAAA VQPAEVTVEV GEDLHMHQVR DREMPEVVEI RRSNCTNHCD LGDTSSYHTK VSTVHIMKKR NGGGSLNNYS SSIPPTPSTS QEDPQFSVPP TANTPTPVCK RSMRWSNLFT SEKGSDPDKE RKAPENHADT

image



Background

Alternative Names	Arf-GAP with GTPase, ANK repeat and PH domain-containing protein 9 (AGAP-9) (Centaurin-gamma-like family member 6)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AGAP9 Protein
Catalog #: PEH10115

