



Product Name: Recombinant Human ZNG1C Protein
Catalog #: PEH10117



Summary

Name	CBWD3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JTY5
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSALKSLAV SQGGELYEEW LELRNGCLCC SVKDNGLRAI ENLMQKKGKF DDILLETGL ADPGAVASMF WVDAELGSDI YLDGIITVD SKYGLKHLTE EKPDGLINEA TRQVALADII LINKTDLVPE EDVKKLRTTI RSINGLGQIL ETQRSRVDS NVLDLHAFDS LSGISLQKKL QHVPGTQPHL DQSIVTITFE VPGNAKEEHL NMFIQNLLWE KNVRNKDNHC MEVIRLKGLV SIKDKSQQVI VQGVHELVDL EETPVSWKDD TERTNRLVLI GRNLDKDILK QLFIATVTET EKQWTTHFKE

image



Background

Alternative Names COBW domain-containing protein 3 (Cobalamin synthase W domain-containing protein 3)

Product Name: Recombinant Human ZNG1C Protein
Catalog #: PEH10117



Note

For Research Use Only , Not for Diagnostic Use