


Summary

Name	C144A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A2RUR9
Host	E.coli
Species	Human
Molecular Mass	69 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TKNMNQNSDS GSTNNYKSLK PKLENLSSLP PDSRTSEVY LHEELQQDMQ KFKNEVNTLE EEFLALKKED VQLHKDVEEE MEKHRSNSTE LSGTLTDGTT VGNDDDGLNQ QIPRKENGEGH DRPADKTSNE KNEVKNQIYP EADFADSMEP SEIASEDCEL SHSVYENFML LIEQLRMEYK DSASLPRIQD TFCLCEHLLK LKNNHCDQLT VKLKQMENMV SVLQNELSET KKTKLQLELQ KIEWEKELYD LRLALKQENE EKRNADMLYN KDSEQLRIKE EECGKVVETK QQLKWNLRRL VKELRTRVNN LDLVVQERN D A QKQLSEEQD ARILQDQILT SKQKELEMAR KKMNSEISHR HQKEKDLFHE DCMLQEEIAL LRLEIDTIKN QNKQKEKKYF EDIEAVKEKN DNLQKIIKLN EETLTETILQ YSGQLNNLTA ENKILNSELE NGKQNQERLE IEMESYRCRL AA AVRDCDQS QTARDLKLDF QRTRQEWVRL HDKMKVDMMSG LQAKNEILSE KLSNAESKIN SLQIQLHNTR DALGRESLIL ERVQRDLSQT QCQKKETE QM

image



Product Name: Recombinant Human CCDC144A Protein
Catalog #: PEH10118



Background

Alternative Names

Coiled-coil domain-containing protein 144A

Note

For Research Use Only , Not for Diagnostic Use