


Product Name: Recombinant Human HERVK_113 Protein
Catalog #: PEH10123



Summary

Name	VP113_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P63121
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VCKAIQ GKQ FEGLVDTGAD VSIILNQWP KNWPKQKAVT GLVGISTASE VYQSTEILHC LGPDNQESTV QPMITSIPLN LWGRDLLQQW GVEITMPAPL YSPTSQKIMT KMGYIPGKGL GKNEDGIKVP VEAQINQERE

image



Background

Alternative Names	Endogenous retrovirus group K member 113 Pro protein (HERV-K113 envelope protein) (HERV-K_19p13.11 provirus ancestral Pro protein) (EC 3.4.23.50) (Protease) (Proteinase) (PR)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HERVK_113 Protein
Catalog #: PEH10123

