


Product Name: Recombinant Human GAB4 Protein
Catalog #: PEH10128



Summary

Name	GAB4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2WGN9
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGHEVQLPPV NRSLKPNQKA NPTPPNLRNN RVINELFKP PVTEPWSGTS HTFDSSSSQH PISTQSITNT DSEDSGERYL FPNPASAFPV SGGTSSSAPP RSTGNIHYAA LDFQPSKPSI GSVTSGKKVD YVQVDLEKTQ ALQKTMHEQM CLRQSSEPPR

image



Background

Alternative Names	GRB2-associated-binding protein 4 (GRB2-associated binder 2-like) (GAB2-like) (GRB2-associated binder 4) (GRB2-associated-binding protein 2-like) (Growth factor receptor bound protein 2-associated protein 4)
--------------------------	---

Note

Product Name: Recombinant Human GAB4 Protein
Catalog #: PEH10128



For Research Use Only , Not for Diagnostic Use