


Product Name: Recombinant Human HEATR5A Protein
Catalog #: PEH10129



Summary

Name	HTR5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86XA9
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CLEDGIRNGS CSSDGLLDLV YADLGTL SRL WLAALQDFAL LTLPSEFASQ LPAEGGAFYT AETSENAKLH YNSWALILH ATALWLTSTG FVVADPDEGA SNLSRPVTPT SMCQGSSSGA TIKSPEDVYT DRFHLILGIS VEFLCSLRSD ATMESITACL HALQALLDVP WPRSKIGSDQ DLGIELLNVL HRVILTRESP SIQLASLEV V RQIICAAQEH VKEKRRSAEV DDGAAEKETL PEFGEGKDTG GLVPGKSLVF ATLELCVCIL VRQLPELNP K LTGSPGVKAT KPQILLEDGS RLVSAALVIL SELPAVCSPE GSISILPTIL

image



Background

Alternative Names HEAT repeat-containing protein 5A

Product Name: Recombinant Human HEATR5A Protein
Catalog #: PEH10129



Note

For Research Use Only , Not for Diagnostic Use