


Product Name: Recombinant Human KLHL6 Protein
Catalog #: PEH10132



Summary

Name	KLHL6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WZ60
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TMGDVVEKSL EGPLAPSTDE PSQKTGDLVE ILNGEKVKFD DAGLSLILQN GLETLRMENA LTDVILCVDI QEFSCHRVVL AAASNYFRAM FCNDLKEKYE KRIIKGVDA ETMHTLLDYT YTSKALITKQ NVQRVLEAAN LFQFLRMVDA CASFLTEALN PENCVGILRL ADTHSLDSLK KQVQSYIIQN FVQILNSEEF LDLPVDTLHH ILKSDDLYVT EEAQVFETVM SWVRHKPSER LCLLPYVLEN VRLPLDPWY FVETVEADPL IRQCPEVFPL LQEARMYHLS GNEIISERTK PRMHEFQSEV

image



Background

Alternative Names Kelch-like protein 6

Product Name: Recombinant Human KLHL6 Protein
Catalog #: PEH10132



Note

For Research Use Only , Not for Diagnostic Use