

Product Name: Recombinant Human LORF2 Protein
Catalog #: PEH10139



Summary

Name	LORF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00370
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TTDPTEIQT IREYYKHLA NKLENLEEMD TFLDTYTLPR LNQEEVESLN RPITGSEIVA IINSLPTKKS

image



Background

Alternative Names	LINE-1 retrotransposable element ORF2 protein (ORF2p) [Includes: Reverse transcriptase (EC 2.7.7.49); Endonuclease (EC 3.1.21.-)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use