


Product Name: Recombinant Human MAGEA9 Protein
Catalog #: PEH10142



Summary

Name	MAGA9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P43362
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPDEDLEAQG EDLGLMGAQE PTGEEEEETS SSDSKEEEVS AAGSSSPQSQ PQGGASSSIS VYYTLWSQFD EGSSSQEEEE PSSVDPAQL EFMFQEALKL KVAELVHLL HKYRVKEPVT KAEMLESVIK NYKRYFPVIF GKASEFMQVI FGTDVKEVDP AGHSYILVTA LGLSCDSMLG DGHSMMPKAAL LIIVLGVILT KDNCAPPEVI WEALSVMGVY VGKEHMFYGE PRKLLTQDWV QENYLEYRQV PGSDPAHYEF LWGSKAHAET SYEKVINYLVLMLNAREPICYPSLYEEVLGE

image



Background

Alternative Names	Melanoma-associated antigen 9 (Cancer/testis antigen 1.9) (CT1.9) (MAGE-9 antigen)
--------------------------	--

Product Name: Recombinant Human MAGEA9 Protein
Catalog #: PEH10142



Note

For Research Use Only , Not for Diagnostic Use