


Product Name: Recombinant Human PLEKHJ1 Protein
Catalog #: PEH10150



Summary

Name	PKHJ1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NW61
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MRYNEKELQA LSRQPAEMAA ELGMRGPKKG SVLKRRLVKL VVNFLFYFRT DEAEPVGALL LERCRVVREE PGTFSSIFIE DPERKYHFEC SSEEQCQEW EALRRASYEF MRRSLIFYRN EIRKVTGKDP LEQFGISEEA RFQLSGLQA

image



Background

Alternative Names	Pleckstrin homology domain-containing family J member 1 (PH domain-containing family J member 1) (Guanine nucleotide-releasing protein x)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use