


Product Name: Recombinant Human PABPC1L2A Protein
Catalog #: PEH10151



Summary

Name	PAP1M_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JQF8
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASLYVGD ^{LH} PEVTEAM ^{LYE} KFSPAG ^{PILS} IRICRD ^{KITR} RSLGYAY ^{VNY} QQPVD ^{AKRAL} ETLNFD ^{VIKG} RPVRIM ^{WSQR} DPSLRK ^{SGVG} NVFIK ^{NLGKT} IDNKAL ^{YNIF} SAFGNIL ^{SCK} VACDEK ^{GPKG} YGFVHF ^{QKQE} SAERAID ^{VMMN} GMFLN ^{YRKIF} VGRFK ^{SHKER} EAERGA ^{WARQ} STSADV ^{KDFE} EDTDEE ^{ATLR}

image



Background

Alternative Names Polyadenylate-binding protein 1-like 2 (RNA-binding motif protein 32) (RNA-binding protein 32)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PABPC1L2A Protein
Catalog #: PEH10151

