



Product Name: Recombinant Human POTEJ Protein
Catalog #: PEH10152



Summary

Name	POTEJ_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P0CG39
Host	E.coli
Species	Human
Molecular Mass	71 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLKLTSEES QRFKGSNSQ PEKMSQEPEI NKDGDREVEE EMKKHESNNV GLENLSNGV TAGNGDDGLI PQRKSRTPEI QQFPDNESEE YHRICELVSD YKEKQMPKYS SENSNPEQDL KLTSEESQR LKGSENGQPE KRSQEPEINK DGDRELENFM AIEEMKKHGS THVGFPENLT NGATAGNGDD GLIPPRKSRT PESQQFPDTE NEEYHSDEQN DTQKQFCEEQ NTGILHDEIL

image



Background

Alternative Names POTE ankyrin domain family member J

Note

Product Name: Recombinant Human POTEJ Protein
Catalog #: PEH10152



For Research Use Only , Not for Diagnostic Use