


Product Name: Recombinant Human PRR20B Protein
Catalog #: PEH10154



Summary

Name	PR20B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P86481
Host	E.coli
Species	Human
Molecular Mass	65 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEEPRPSKRL RSMAPNQASG GPPPEPGCCV ADPEGSVEAD GPAQPAQPAK PIAYVKPFRR QPPARPESPP PAERGRRRGG SRRPGRGRGR RAGPRGDAGQ RQGAEGLMAP DVHIQLDHHG EPGHQGEPEI TETAAFLSE TGPPPQTVQE GPGPDVAQPE LGFQEPPAAP GPQAVDWQPV LTLYPCIGFR ALGDSAVLQV IQTPQGTYVQ GVPVFLTDIA Y

image



Background

Alternative Names Proline-rich protein 20B

Note

Product Name: Recombinant Human PRR20B Protein
Catalog #: PEH10154



For Research Use Only , Not for Diagnostic Use