


Product Name: Recombinant Human FAM47B Protein
Catalog #: PEH10159



Summary

Name	FA47B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NA70
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DAWARCEARE KTTEVPTESG KYPCGESCPR PPETPVSRLR PQLPKTPVSS RRPEPPKTRV SSLRPEPPKT RVSSLHPEPP ETRASHLRVD PPETGVSHLC PEPPKTLVSS VHPEPPDTGA SHLCPEPPET RVSHLHPEPP ETGVSHLRPE PSKTQVSSLC PEPPEAGVSH LCLEPPNTHR VSSFLLQVLK LDSEKKLEDA RARCEGQEMT

image



Background

Alternative Names Protein FAM47B

Note

Product Name: Recombinant Human FAM47B Protein
Catalog #: PEH10159



For Research Use Only , Not for Diagnostic Use