


Product Name: Recombinant Human CROCCP3 Protein
Catalog #: PEH10167



Summary

Name	CROL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVE0
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSESERRALE EQLQRLDKT DSTMQAHEDA QREVQRLRSA KELLRREKSN LAHSLQVAQQ QAKELRQERK KLQAAQEELR RQRYWLGEELQ EDAVQDGVRV RRELERSHRQ LEQLEGKRSV LAKELVEVRE ALSRATLQRD MLQAEKAEVA

image



Background

Alternative Names	Putative ciliary rootlet coiled-coil protein-like 2 protein (Ciliary rootlet coiled-coil protein pseudogene 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use