


Product Name: Recombinant Human FAM90A26 Protein
Catalog #: PEH10180



Summary

Name	F90AP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	D6RGX4
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEKEERPRQQ DPQRKALLHI FSGKPPEKPL PNRKGSTESS VYLRVASGPM PVHTTSKRPR VDPVLADRSA TEMSDRGSAL ASLSPLRKAS LSSSSSLGPK ERQTGAAADI PQPAVRHQGP EPLLVKPTH SSPEGGCREV PQAASKTHGL LQAISPQAQD KRPAVTSQPC

image



Background

Alternative Names Putative protein FAM90A26

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FAM90A26 Protein
Catalog #: PEH10180

