


Product Name: Recombinant Human FAM90A2P Protein
Catalog #: PEH10181



Summary

Name	F90A2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q658T7
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NLKPWKQVE ANGPLNKDK GEKEERPRQQ DPQRKALLHI FSGKPPEKLL PNRKGSTESS VFLRVASRPM PVHTTSKRPC VDEPADDRSA TEMSGRGSVL ASLSPLRKAS LRSSSLGPK ERQTGAAADI PQPAVRHQGP EPLLVKPTH SSREGGCREV PQAASKTHGL LQAVRPQAQD KRPAVTQPGP PAATHSLGLG

image



Background

Alternative Names Putative protein FAM90A2P

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FAM90A2P Protein
Catalog #: PEH10181

