


Product Name: Recombinant Human LINC00315 Protein
Catalog #: PEH10186



Summary

Name	CU093_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P59091
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDSLTERCRP QPLAALGPAD AQLHGPAAG GGVGPRVMLS SLCREHPWQG RCPPIPAAGK EVRQLRHLHR TRPPDTHQDM ARGPLPHTH GLPARPSHVR HQTEPVQVAS RWHVPSDGLC RPCHVHQLLH PPLVPALGA

image



Background

Alternative Names Putative uncharacterized protein encoded by LINC00315

Note

For Research Use Only , Not for Diagnostic Use