


Product Name: Recombinant Human MIR1915HG Protein
Catalog #: PEH10189



Summary

Name	CSC10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T4H9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TAERRRGWRC RGVPTPSGDR GPGAARPAAR GGAGETTQQ RPLQVPGASA AAARTLLRW HHRVPSPRAT RSPGSIRRTS PCSGGPDRPE

image



Background

Alternative Names	Putative uncharacterized protein encoded by MIR1915-HG (Cancer susceptibility 10 protein) (Cancer susceptibility candidate gene 10 protein) (MIR1915 host gene protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use